


# Materials Science: Solar Conversion Efficiency

## Correlation between photoluminescence and carrier transport and a simple *in situ* passivation method for high-bandgap hybrid perovskites.

R.J. Stoddard, F.T. Eickemeyer, H.W. Hillhouse et. al., J Phys Chem Lett (2017) **8**: 3289–3298.

- Hybrid perovskites (HPs) have outstanding optoelectronic characteristics for use in photovoltaic devices.
- Measurements of photoconductivity and photoluminescence quantum yield (PLQY) are used for assessment of HP performance.
- Lumencor's MAGMA Light Engine delivers unprecedented performance by integrating 21 individually selectable, bright, stable, solid-state light sources across the UV-VIS-IR spectrum. Spectral content may be tuned by using several independent color band outputs in series or simultaneously, all under the control of an onboard microprocessor.
- Sources are well aligned, provide uniform illumination and color temperature control. Outputs are managed for robust, reproducible operation over time within one Light Engine and between multiple Light Engines.
- Hyper-spectral imaging, solar simulation and various metrology instrumentation applications are supported.



## MAGMA Light Engine for Solar Energy Conversion

Experimental setup shown schematically for simultaneous measurement of the PLQY and mean carrier diffusion length.

Calibrated LED excitation, quantitative wide-field photoluminescence and four-point photoconductivity is used.

The blue circle indicates the illuminated area, (excitation with a blue light source from a Lumencor MAGMA Light Engine.)

Dashed lines indicate scribe lines to avoid spastic current pathways.

The current is sourced along gold (Au) contact pads 1 and 2, and the voltage is measure with inner probes (between 3 and 4 or 5 and 6).