AURA II Light Engine[®]

Instruction Manual

Regulatory Models

Lumencor utilizes regulatory model names for all certified and CE marked products. The regulatory model names are traceable to all regulatory documentation, third party reports and certifications.

“Regulatory Model: Aura II” is used as a representative model for all certified and CE marked Aura II Products.

Emissions

This equipment has been tested and found to comply with the limits of EMC directive 2014/30/EU. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Safety Certifications

TUV SUD America, CB Certification (IEC 61010-1:2010)

TUV SUD America, NRTLus Certification (UL 61010-1:2012-05)

TUV SUD America, cNRTL Certification (CAN/CSA-C22.2 No. 61010-1:2012)

TUV SUD America, EN Certification (EN 61010-1:2010)

CE Marking

Low Voltage Directive (2014/35/EU)

EMC Directive (2014/30/EU)

RoHS Directive (2011/65/EU+2015/863/EU)

REACH Regulation (EC) No. (1907/2006/EC)

EU Declarations of Conformity can be found at <https://lumencor.com/company/compliance>

For EU customers discarding end-of-life Lumencor electrical and electronic equipment: Please submit an RMA request with “Recycle product under WEEE” in the Description of Issues field.

For disposal in countries outside of the European Union: This symbol is only valid in the European Union (EU). If you wish to discard this product, please contact your local authorities or dealer and ask for the correct method of disposal.

Lumencor Light Engines as supplied, and as represented in this manual, meet safety and regulatory requirements For Research Use Only. If the light engine is incorporated into an instrument or system for a specific end-use application, it is the responsibility of the system integrator to verify that the light engine, and the system into which it is incorporated, meet all safety and regulatory requirements of that end-use application.

Lumencor, Inc. | 14940 NW Greenbrier Parkway | Beaverton, OR 97006 USA | 503.213.4269 | www.lumencor.com

Document Number 57-10007

Revision B

092121

Table of Contents

1. Introduction
2. Precautions and Warnings
3. Installation and Operating Instructions
4. Spectral Output
5. Product Specifications
6. Routine Maintenance and Trouble Shooting
7. Customer Support
8. Warranty

1. Introduction

Lumencor AURA II light engines are designed for laboratory use by bioanalytical researchers and/or developers of bio-optical instrumentation. The AURA II provides 2, 3, 4, or 5 individually controllable solid state light sources that have been optimized to produce spectrally discrete outputs. This manual covers all AURA II models identified by the designation **Regulatory Model: AURA II** on the identification label, located on the back panel of the light engine (see Figure 1).

The light sources within the AURA II are controlled by a serial command interface connecting to the Lumencor light engine control pod accessory (83-10007) or a computer running a Lumencor supplied GUI or a third party microscopy software application. Alternatively, the light sources may be enabled by TTL inputs from a trigger device such as a camera. The user can enable or disable each source independently (serial or TTL) and change the intensity of each source independently (serial only). The only manual control is a power switch on the front panel to turn on/off the electrical power supply to the unit. A green indicator above the power switch is lit when the power supply is connected to the AURA II and the power switch is in the on position.


Figure 1. Aura II identification label.

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A few simple practices will ensure trouble-free operation for the life of the light engine.

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Safety Instructions:

Please read and follow all safety instructions provided **BEFORE** using your new AURA II. Failure to comply with the safety instructions may result in fire, electrical shock, or personal injury and may damage or impair protection provided by equipment. Please save all safety instructions.

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Warning: Statements identify conditions or practices that could result in personal injury.

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Caution: Statements identify conditions or practices that could result in damage to your equipment.

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Safety Items {Mesures de sécurité}:

Warning: ONLY use the power supply provided by Lumencor. The Lumencor-supplied 24 VDC, 9.2 A external power supply is required for use with the AURA II light engine. It is imperative that the DC power supply has output over-current protection, as the power input of the AURA II is not fused. The DC power supply must have the AC power cord connected to a receptacle with a protective safety (earth) ground terminal.

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Warning: DO NOT stare into the output of the light engine. The brightness of this light source is higher than most commercial lighting fixtures and is intended to couple directly into a microscope or other bioanalytical instrument.

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Warning: DO NOT turn on the light without the output end of the light guide safely directed into an enclosed optical path. DO NOT point the light output directly onto any flammable or burn-susceptible material. This includes all animal or vegetable tissues, plastics, fabrics, paper and liquids.

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Warning: Possibly hazardous optical radiation emitted from this product. Do not look at operating lamp. Eye injury may result.

Warning: Infrared (IR) emitted from this product. Do not look at operating lamp.

Warning: UV emitted from this product. Avoid eye and skin exposure to unshielded product.

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Caution: DO NOT open the unit. There are no serviceable parts inside and opening the light engine enclosure will void the manufacturer's warranty.

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Caution: DO NOT connect a video cable to the TTL input enable port. Although the connector might look compatible, this input is not intended to be driven by a video signal.

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Caution: DO NOT set liquids on the light engine. Spilled liquids may damage your light engine.

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Caution: DO NOT drop the light engine. It contains glass optical components that could be damaged or misaligned by the shock produced by a drop onto a hard surface.

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DISCLAIMER: Lumencor shall not be liable for injury to the user or damage to the product resulting from the AURA II being used in a way for which it was not intended and in complete disregard for any posted safety precautions and warnings.

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3. Installation and Operating Instructions

All AURA II family light engines ship with the following list of standard components:

1. AURA II light engine, configured with three, four or five solid-state light sources (color channels). The output of each source is refined by a bandpass filter. Bandpass filters are installed to order according to customer specifications and are not user-exchangeable. The filter configuration for each AURA II light engine is shown on the certificate of conformance issued to the customer when the light engine is shipped by Lumencor (Figure 2).
2. A 24V/9.2A (220 W) DC power supply (Lumencor part no. 27-10019)
3. A 6ft AC power cord for the power supply (for North American customers, Lumencor part no. 29-10002, for UK customers, Lumencor part no. 29-10004 and for European customers, Lumencor part no. 29-10005).
4. A USB A-to-RS232 cable (Lumencor part no. 29-10011) **or** USB B (M) to USB A (M) cable (Lumencor part no. 29-10058)


Figure 2. Example certificate of conformance, AURA II Light Engine.


Figure 3. Aura II Light Engine front panel.

NOTE: Any end-product/system incorporating or coupled to a Lumencor Light Engine shall be fully evaluated to verify all applicable safety and regulatory compliance requirements prior to use.

AURA II light engines incorporate fan-assisted air cooling. When positioning the AURA II for use, place the unit on a hard surface and avoid blocking or restricting airflow at the air intake (front panel) and exhaust (rear panel) ports. The AURA II must be operated in a non-condensing environment (dew point $<10^{\circ}\text{C}$ with controlled ambient temperature $<30^{\circ}\text{C}$). Restricting the airflow will cause the unit to operate at elevated temperatures and will result in decreased operating lifespan and/or premature failure.

Warning: DO NOT turn on the light without the output safely directed into an enclosed optical path.

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AURA II light engines may be configured for light delivery via a 3 mm liquid light guide or an SMA-terminated optical fiber. The output couplers are built into the light engine front panel and are not field-exchangeable (Figure 3). The 3 mm liquid light guide output coupler incorporates an electrical power interlock switch in the back of the receptacle. The liquid light guide must be fully inserted in the receptacle to activate the interlock. Electrical power supply to the light engine cannot be turned on via front panel power switch until the interlock is activated. After fully inserting the liquid light guide, lock it in position using the set screw on the side of the receptacle to prevent inadvertent disconnection of the light guide during use.

The AURA II can be controlled via serial commands delivered through the RS-232 port from the Lumencor light engine control pod accessory (83-10007), a Windows GUI (free download from Lumencor's [Control Software site](#)) or via customer-supplied instrumentation control software. Serial commands provide controls for on/off switching and intensity adjustment (0 to 100% in 1% increments) of individual color channels. Refer to Figure 3 for the location of the various connectors. TTL signals delivered to the DB15 connector (Figure 3) provide high speed (10 μs) on/off switching of individual color channels.

3.1 GUI Installation and Operation

1. Operation and installation of the AURA II GUI requires a computer running the Windows operating system with a free USB port.
2. Download the zip file for the AURA II GUI from Lumencor's [Control Software site](#).
3. Unzip the file and run setup.exe to install the AURA II GUI.
4. Connect the USB-to-RS-232 cable from the computer to the DB9 RS-232 port on the AURA II (Figure 3) [1].
5. Successful installation is indicated by the appearance of "USB Serial Port (COM #)" under the "Ports (COM & LPT)" tab in the Windows Device Manager. If the virtual COM port (VCP) is not registered by the operating system, download and install the VCP driver from Lumencor's [Control Software site](#).
6. Connect the DC power supply to the AURA II.
7. Check that the liquid light guide is fully inserted and locked in the front panel receptacle (Figure 3).
8. Toggle the power switch on the front panel to the ON position. The green LED above the switch (Figure 3) should light.
9. Run the GUI by going to the Program Menu and selecting LLE Controller.
10. In the COM pulldown menu (GUI window, upper left), select COM # assigned to USB-Serial port.
11. Press the "INIT" button in the GUI
12. The computer should now have control of the AURA II. Graduated sliders control the source output intensity for each color channel. ON/OFF buttons for each color channel are located below the respective slider controls [2, 3].

Notes

- [1] On some AURA II models, the port is USB B and the connecting cable is USB B (M) to USB A (M).
- [2] The identity of the source outputs corresponding to control channels 1–5 differs among AURA II light engine models. If you are unsure of the correspondence, please contact our Technical Support Group (see Section 7). Please provide notification of the serial number of your AURA II light engine with your inquiry.
- [3] Two luminescent light pipe sources (e.g. Green and Yellow, Figure 5A) cannot be ON simultaneously.

Connector Pin Definitions

DB15HD connector (TTL)		
Pins	Definition	DC Characteristics
1, 2, 3, 12, 15	Channel 5, 3, 2, 1, 4 enables respectively.	VCC = 5.0V V _{ilow} (max) = 1.5V, V _{ihigh} (min) = 3.3V, I _{ilow} = 0.5mA, I _{ihigh} = 1.0μA
6, 7, 8, 10	Gnd	
5	Dosimeter	
4, 9, 11,13,14	N/C	

DB9 connector (RS-232)		
Pins	Definition	DC Characteristics
1, 2, 3, 4, 6, 7, 8	DCD, RXD, TXD, DTR, D3R, RTS, CTS	VCC = 5.0V V _{ilow} (max) = 0.8V, V _{ihigh} (min) = 2.4V, I _{ilow} = 0.5mA, I _{ihigh} = 1.0μA
5	Gnd	
9	N/C	

3.2 Light Engine Control Pod Interface

1. Connect the USB A port of the light engine control pod accessory (83-10007) to the DB9 port on the AURA II (Figure 3) using the USB-to-RS-232 cable (29-10011) [1]. Connect the USB B port of the control pod to a USB A port on a host computer. For applications not requiring a host computer, connect a DC power supply (5V/ 500mA) to the USB B port. Press and hold the **right** button on the pod until a menu of light engines appears.


Figure 4. USB connectors on rear of the light engine controller pod

Turn the rotary dial to select “AURA” from the menu. Press the right button again to return to the main (0–100 analog intensity) display screen.

2. Press the **left** button to select the desired color channel. Successive presses will cycle through the available color channels.
3. Press the **right** button to turn the selected light source on. Press the **right** button again to turn the selected light source off. Adjust the output intensity using the rotary dial [2,3].
4. Press and hold the **left** button to obtain digital output intensity settings for each color channel [4]. Press the **right** button to return to the main display screen.

Notes

- [1] On some AURA II models, the port is USB B and the connecting cable is USB B (M) to USB A (M).
- [2] Output intensity can be set from 0–100% in 1% increments; however operation in the 0–5% range is not recommended.
- [3] Dialing intensity to zero automatically issues an OFF command to the light engine. Press the **right** button to turn the light output on again.
- [4] The current intensity settings are internally stored. When the pod is powered down, the settings are retained and will be restored at the next restart.

3.3 TTL Interface

The TTL Interface provides users with a faster method of switching color channel outputs on and off. Individual TTL lines are provided for each color channel. These can be conveniently addressed using an accessory BNC breakout cable (Lumencor part no. 29-10015 or 29-10080) connected to the front panel DB15 port (Figure 3). As a safeguard against unintended light output when the inputs are initially connected, the TTL port is disabled and must be enabled by a serial command. To enable the TTL port, establish a serial connection as described in Section 3.1 and click the TTL ON (SFT) button in the GUI. A description of the serial command string for this operation is available on request from Lumencor Technical Support (Section 7). Note that TTL and serial on/off commands have a logical OR relationship. Therefore on/off controls in the GUI or other serial control software must be set in the “OFF” state when using TTL control. TTL polarity (ACTIVE = HIGH or ACTIVE = LOW) is factory set according to order specifications. TTL polarity can be reversed using controls provided in the GUI.

4. Spectral Output

The spectral output of AURA II light engines depends on the configuration of solid-state light sources installed and the bandpass filters associated with each source. Two typical examples are shown in Figure 5. The light sources installed and the bandpass filters associated with each


Figure 5. Output spectra for 4-channel AURA light engines

are shown on the Certificate of Conformance included with the shipping documents e-mailed by Lumencor to the customer. The certificate of conformance also shows full (100%) power outputs for each source measured at the terminus of the liquid light guide (LLG) or optical fiber (Figure 2).

5. Product Specifications

AURA II light engines must be operated and stored within the environmental conditions specified.

Specification	Detail
Temperature	
Operating	32 to 95° F (0 to 35° C)
Non-operating	-4 to 158° F (-20 to 70° C)
Humidity	
Operating and non-operating	0 to 80% relative humidity, non-condensing
Altitude	
Operating	0 to 10,000 feet (3,048 meters)
Non-operating	0 to 45,000 feet (13,176 meters)
Dimensions	
Size (W x L x H)	12.5 cm x 26.3 cm x 16.3 cm (4.9 in x 10.4 in x 6.4 in)
Weight	3.6 kg (8 lb)
Input Power Requirements	24V DC / 7.9A
Warm-up Period	1 s
Ingress Protection	IP Rating of X0
Sound Level	Sound level at 1 meter < 65db(A)
Control Interfaces	Serial (USB to RS-232), TTL
Warranty	36 months parts and labor

Performance specifications for individual light engines are listed on the certificate of conformance included with the shipping documents e-mailed to the customer (see example shown in Figure 2). It is important to retain the certificate of conformance for reference. In the event that the light engine is sold, the certificate of conformance should be transferred to the new owner. Certificates of conformance are also recorded in Lumencor's database and copies can be requested by e-mail to techsupport@lumencor.com. The request message must include the 4- or 5-digit serial number of the light engine.

6. Routine Maintenance and Trouble Shooting

No routine maintenance is required. There are no user-replaceable components or sub-assemblies in the AURA II light engine. Opening the light engine enclosure will void the manufacturer's warranty. In the event that the light engine fails to perform in accordance with the specifications listed on the certificate of conformance, please contact Lumencor Technical Support for assistance, as directed in Section 7.

7. Customer Support

For technical support of the AURA II light engine, please contact Lumencor by phone at 503.213.4269 or via e-mail at techsupport@lumencor.com. Please be prepared to provide the 4- or 5-digit serial number of the light engine (see Figure 1), a description of the problems encountered and information on the usage context (e.g. what microscope and what control software is being used). This information will help to determine whether the problems can be resolved in situ by adjustments to the system configuration, or whether a fault has developed in the light engine that requires its return to Lumencor's facility in Beaverton, Oregon for evaluation and, if necessary, repair. Any light engine return to Lumencor for service or repair requires a material authorization (RMA) number. To request an RMA number, fill out and submit the [online request form](#). It is the customer's responsibility to properly package and safely ship products to Lumencor. Instructions for shipping will be provided in the e-mail giving notification of the RMA number.

8. Warranty

The AURA II family of light engines is backed by a 36 month warranty. An extended warranty option is available; extended warranties must be purchased in advance of receipt of goods or within 90 days of the date of original shipment from Lumencor. Accessories including (but not limited to) liquid light guides, optical fibers, collimators, cables and control consoles are not covered by the warranties attached to light engines. Please fill out and submit the [online warranty registration form](#). This will facilitate provision of warranty service should it be required.